

REVEAL[®] Q+ FOR OCHRATOXIN

DOWNLOAD AND READ KIT INSTRUCTIONS COMPLETELY BEFORE PERFORMING TEST.



Materials Provided:

25 Reveal[®] Q+ for Ochratoxin test strips
25 Red sample dilution cups
25 Sample cups
01 Bottle of sample diluent

Product Number: 8685

Threshold: 2–20 ppb

Testing time: 9 minutes

Call 800.234.5333 to order or visit NEOGEN.com

Procedure



1 Prepare by entering QR code into the AccuScan[®] Gold reader.



2 Obtain a representative sample. Grind and weigh out a 10 g sample.



3 Add 40 mL of methanol to the sample.



4 Shake vigorously for 3 minutes, or blend for 1 minute.



5 Settle, then filter.




6 Transfer 200 μ L to the red dilution cup.



7 Add 100 μ L sample extract to the red dilution cup and mix up and down 5 times.




8 Transfer 100 μ L to sample cup.



9 Place a new Reveal for Ochratoxin strip into the sample cup. Set a timer for 9 minutes.



10 Remove promptly at 9 minutes and interpret results using the AccuScan Gold reader.



Results shown on the AccuScan Gold reader.



Reveal[®] Q+ for Ochratoxin

Aqueous Extraction

Product Number: 8685

Ordering Information

8685 Reveal Q+
for Ochratoxin
AccuScan[®] Gold



Test with Confidence

NEOGEN[®] introduces Reveal[®] Q+ for Ochratoxin, an easy-to-use, rapid, quantitative lateral flow test that provides unparalleled accuracy in only nine minutes.

Reveal Q+ for Ochratoxin is innovative technology which yields accurate and reproducible quantitative results.

- **Simple** water based extraction
- **Fast** time to result
- **Easy** to use

Materials Recommended, Not Provided

NEOGEN #	Item Description
9428, 9428B	Sample collection cups with lids 125 mL
9421, 9421B	Sample collection tubes with caps
9420, 9429, 9519	Filter syringes, Whatman #4 filter paper or equivalent
9330	Centrifuge, mini
9172	Micro centrifuge tubes
9448, 9447	Dispensing pump or graduated cylinder
9401, 9453	Agri-grind grinder or equivalent
9427	Scale capable of weighing 5–50 g ± 0.1 g
9426, 9452	Timer
9475	Reveal sample cup rack
9860, 9272, 9278	Pipettor, 100 µL
9464, 9487	Pipette tips 1–200 µL
9488	Pipettor, 200 µL
—	Distilled or deionized water

